


CRIC

**CRICOTHYROIDOTOMY
MODULE | MP0002**



Turn the V2 into a cricothyroidotomy trainer for both needle and surgical techniques.

The module includes palpable landmarks, a mobile trachea and replaceable skin and tissue inserts.

Change the inserts to vary the difficulty and reset the trainer quickly between attempts.

KEY FEATURES



***MOBILE TRACHEA**

The trachea moves beneath the skin and must be located, stabilised and controlled during the procedure.



***ADJUSTABLE DIFFICULTY**

Use standard or thick tissue inserts to change how easy the airway and cricoid membrane are to palpate.



***SIMULATED BLEEDING**

Add the blood-filled insert to introduce simulated bleeding and make the procedure more challenging.

SET UP AND RESET



FIT THE MODULE

Swap the standard airway and face for the Cricothyroidotomy Module.



CHANGE THE TISSUES

Swap the skin and tissue inserts between attempts.

SUPPORTS THE FOLLOWING LEARNING OBJECTIVES:

ANATOMY

Locate the cricothyroid membrane using palpable landmarks.

STABILISATION

Stabilise a mobile trachea during the procedure.

TECHNIQUE


Practise needle and surgical cricothyroidotomy.

DIFFICULTY

Adapt the procedure when increased tissue depth makes landmarks difficult to palpate.

KIT INCLUDES:

- * Cricothyroidotomy face skin
- * Cricothyroidotomy airway
- * 6 replaceable skin inserts
- * 2 standard tissue inserts
- * 2 obese tissue inserts
- * 1 morbidly obese tissue insert
- * 1 blood-filled insert
- * Repair tape
- * Instructions for use



WATCH THE
CRICOTHYROIDOTOMY →
MODULE IN USE



CONTACT US VIA
[DECENTSIMULATORS.COM](https://www.decentsimulators.com)